

## Chromatin Immunoprecipitation (ChIP)

1. Add **3 ug** antibody and 15 ul of Dynabeads to cross-linked, fragmented chromatin extract.
2. Cap tubes and place on rotator at 4°C 2 hrs to O/N

## ChIP Wash

1. Place sample plate on the magnet for ~1 min.
2. While still on the magnet, aspirate and discard supernate.
3. Wash beads with **NaCl 250 wash buffer**.
  - a. Remove tubes from magnet, and add ~150 ul **NaCl 250** (4°C) to each well.
  - b. Cap the tubes and mix by inverting ten times.
  - c. Centrifuge, 1 sec.
  - d. Place sample plate on the magnet for ~1 min.
  - e. Remove caps. *Keep them in order, and in same orientation.*
  - f. While still on the magnet, aspirate and discard extract.
    - i. *Take care to stay away from beads. No need to change tips. Quickly go to next step.*
4. Wash beads with **LiCl 250 wash buffer**.
  - a. Repeat process described above.
5. Wash beads with 10 mM Tris-HCl
6. Remove from magnet, and immediately add **Polishing mix** (see below)