

ChIP DNA TruSeq Adapter Ligation

1. Fill out **Master Ligation mix** table. **X** = (1.1) x number of ChIP samples

<i>Master Ligation mix (make on ice)</i>	<i>1x (ul)</i>	<i>[Final]</i>
ddH ₂ O	15	
(Enzymatic) 10x T4 DNA ligase buffer	2	1x
(Enzymatic) 600 U/ul T4 DNA Ligase	1	600 U
2mM TruSeq Fork Adapter	2	200 uM
<i>Total reaction volume</i>	(20)	

2. Make **Master ligation mix** in 1.5 ml tube (or 15 ml tube if needed). Vortex to mix.
3. Transfer **20 ul** of **Master ligation mix** to each ChIP samples.
4. In the thermomixer, incubate samples for **2 hrs @ 25°C**.
5. Place tubes against magnet for 1 min., remove caps keeping them in order, draw off supernate.