

**ChIP Elution Buffer (50 mL)**

|                        |       |         |
|------------------------|-------|---------|
| 25 mM Tris-HCl, pH 7.5 | 1 M   | 1.25 ml |
| 2 mM EDTA, pH 8.0      | 0.5 M | 0.2 ml  |
| 200 mM NaCl            | 5 M   | 2 ml    |
| 0.5% SDS               | 20%   | 1.25 ml |

Bring final volume to 50 mL with ddH<sub>2</sub>O, filter and store at RT.