

## ChIP DNA A-tailing

1. Fill out **A-tailing mix** table. **X** = (1.1) x number of ChIP samples

<b>A-tailing mix</b>	<b>1x (ul)</b>	<b>[Final]</b>
ddH <sub>2</sub> O	25	
10x NEB2 buffer	3	1x
<b>3mM dATP</b>	1	100 uM
5 U/ul Klenow 3'-5' exo minus	1	5 U
<i>Total reaction volume</i>	(30)	

2. Make **A-tailing mix** in 1.5 ml tube. Vortex to mix.
3. Transfer **30 ul** of **A-tailing mix** to each sample, and mix gently up and down upon transfer.
4. In the thermomixer, incubate samples for **30 min @ 37°C**.
5. Place tubes against magnet for at least 1 min., remove caps keeping them in order, draw off supernate.