

<b><u>Modified Farnham Lysis Buffer (500 mL)</u></b>	<b><u>Stock</u></b>	<b><u>Volume</u></b>
20mM Tris pH 8.0	1M	10 mL
85mM KCl	2M	21.25 mL
0.5% NP-40 (IGEPAL)	10%	25 mL
0.5% Triton X-100	25%	10 mL

Bring final volume to 500 mL with ddH<sub>2</sub>O. Filter sterilize, store at 4°C.

Source: <https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0100099#pone-0100099-t001>